

Positive Control Cell Lysate

BSA-Nitrotyrosine

Order No.: 0068-2000/NOY-BSA

Lot: 0068



www.nanotools.de

orders & support:

email: info@nanotools.de

phone: +49-7641-455 670

fax: +49-7641-455 671

03/221013F

Formulation Nitrotyrosine - modified BSA
4 mg/ml in PBS/0,09% Na - azide.

Reconstitution Aliquote and freeze in liquid nitrogen; NOY-BSA can be stored frozen at -20°C up to 1 year. Thaw aliquots at 37°C.

Avoid repeated freeze/thaw cycles.

Application BSA-Nitrotyrosine is recommended as positive control for immunoblot applications. Use 50 ng per lane for HRPO/ECL detection.