

Positive Control Cell Lysate

Neuro 2A untreated

Order No.: **0911/liquid/Neuro2A/Co**

Lot: **0911**



www.nanotools.de

orders & support:

email: info@nanotools.de

phone: +49-7641-455 670

fax: +49-7641-455 671

02/260207F

Formulation

Cell lysate from serum starved Neuro 2A cells diluted with SDS-PAGE sample buffer.

Reconstitution

Aliquote and store frozen.

Avoid repeated freeze/thaw cycles.

Application

The positive control cell lysate is recommended for immunoblot applications. 20µl of positive control cell lysate correspond to ca. 20.000 cells. Heat control cell lysate 5min at 90°C before use.

Use 20µl / lane (mini gel) for HRPO/ECL detection of the target proteins.