

Positive Control Cell Lysate

SKOV-3 untreated

Order No.: 0991/SKOV3/Co

Lot: 0991



www.nanotools.de

orders & support:

email: info@nanotools.de

phone: +49-7641-455 670

fax: +49-7641-455 671

02/260207F

Formulation

Lyophilized cell lysate from serum starved SKOV-3 cells.

Reconstitution

Reconstitute by addition of 200 μ l H₂O. After complete solubilization add 200 μ l 2x SDS-PAGE sample buffer, mix and incubate at 90°C for 5 min.

Aliquote and store frozen. Avoid repeated freeze/thaw cycles.

Application

The positive control cell lysate is recommended for immunoblot applications. 20 μ l of positive control cell lysate correspond to ca. 20.000 cells.

Use 20 μ l / lane (mini gel) for HRPO/ECL detection of the target proteins.

Please note:

The lyophilized cell lysates contain SDS and **are not recommended** for applications with native proteins such as immunoprecipitation.