

## Positive Control Cell Lysate

### SH-SY5Y pervanadate treated

Order No.: 1002/SH-SY5Y/VH

Lot: 1002/



[www.nanotools.de](http://www.nanotools.de)

orders & support:

email: [info@nanotools.de](mailto:info@nanotools.de)

phone: +49-7641-455 670

fax: +49-7641-455 671

02/260207F

### Formulation

Lyophilized cell lysate from SH-SY5Y cells.  
Serum starved cells were treated for 15min with pervanadate.

### Reconstitution

Reconstitute by addition of 200  $\mu$ l H<sub>2</sub>O. After complete solubilization add 200  $\mu$ l 2x SDS-PAGE sample buffer, mix and incubate at 90°C for 5 min.

**Aliquote and store frozen. Avoid repeated freeze/thaw cycles.**

### Application

The positive control cell lysate is recommended for immunoblot applications. 20 $\mu$ l of positive control cell lysate correspond to ca. 20.000 cells.

Use 20 $\mu$ l / lane (mini gel) for HRPO/ECL detection of the target proteins.

Please note:

The lyophilized cell lysates contain SDS and **are not recommended** for applications with native proteins such as immunoprecipitation.