

Immunoblotting Reagents

Blot Storage Buffer BSB

Order No.: 3030-0500/BSB
Lot: 3030B
Size (ml) 500



www.nanotools.de

orders & support:

email: info@nanotools.de
phone: +49-7641-455 670
fax: +49-7641-455 671

03/080507F

Formulation

The Blot storage buffer BSB is a PBS / Casein buffer containing detergents. BSB contains phosphatase inhibitors and components that protect membrane bound proteins during long term storage at -20°C.

Storage

For long term storage, freeze buffer (-20°C) upon arrival. Thawed buffer may be stored at 4°C up to 1 week.

Application

Use 1x blot incubation buffer to dilute the primary and the secondary antibody to the appropriate concentration.