

## Immunoblotting Reagents

### Blot Incubation Buffer CPPT

Order No.: 3031-0500/CPPT  
Lot: 3031B  
Size (ml) 500



[www.nanotools.de](http://www.nanotools.de)

**orders & support:**

email: [info@nanotools.de](mailto:info@nanotools.de)  
phone: +49-7641-455 670  
fax: +49-7641-455 671

03/080507F

### Formulation

The Blot incubation buffer CPPT is a PBS / Casein buffer containing detergents. The buffer does not contain preservatives.

### Storage

For long term storage, freeze buffer (-20°C) upon arrival.  
Thawed buffer may be stored at 4°C up to 1 week.

### Application

Use 1x blot incubation buffer to dilute the primary and the secondary antibody to the appropriate concentration.